

Characterisation of Bi-refrangent coatings and substrates using the Aquila nkd Spectrophotometer, Optical Coating Measurement System.



This application note illustrates the use of the nkd- series optical coating measurement system for characterising transparent bi-refrangent films or substrates with ease.

Bi-refrangent or optically anisotropic materials have optical properties which are not the same in all directions within a given sample. Some natural crystal materials exhibit birefringence and other materials such as polymers can have birefringence induced intentionally or not by the manufacturing process. Stretched polymer films can have different refractive indices for different orientations of the sample. The difference in refractive index for s- and p-polarised incident light, is used as a measure of the birefringence of the material.

Measuring transparent bi-refrangent films and substrates poses an impossible task for old-fashioned spectrophotometers and a difficult one for ellipsometers. For the old-type spectrophotometers, the difficulty arises from the fact that the incident light polarisation is unknown and undefined. Without defining the incident polarisation, no sensible analysis of the resulting T and R spectrum can be performed. The Aquila nkd spectrophotometer overcomes this because the incident light polarisation is well known and can be selected for each measurement.

In ellipsometry the incident polarisation is known but other problems arise from the fact that the samples are transparent. Ellipsometric analysis cannot account for the multitude of front and back surface reflections generated by transparent substrates. To overcome this the sample needs to be roughened on one side making this a destructive measurement technique. Even with this roughening or blacking out treatment, it is never an ideal measurement and the valuable transmission data is lost. With the nkd spectrophotometer, transparent substrates are an advantage as the additional spectroscopic information provided by the transmission data allows for accurate determination of the sample absorption. No sample preparation is required and the user simply places the film or substrate on the sample platform. The incident polarisation can be selected using the Pro-Optix™ software, which handles all aspects of the measurement and analysis process.

Requirement for measuring these anisotropic materials is increasing with the processing of thin glasses and polymer films for display applications. LCD image quality for example can suffer if there is spatial variation in light polarization over the screen, producing dark streaks on the viewing area. The birefringent effect can also be very beneficial. Liquid crystals are a twisted birefringent material used in many laptop displays. The example here examines one type of polymer film which has been processed intentionally to induce bi-refringence.

In Figure 1, a sample of Polycarbonate stretched in one direction has been measured on the nkd-7000. Transmittance (T) and reflectance (R) measurements for s- polarised incident light are shown. A note of the orientation of the sample for this measurement was made and is indicated as 0° sample orientation in the results. Analysis of this spectra using the powerful Pro-Optix™ yielded the dispersion shown in Figure 2. A Cauchy model, just one of the dispersion algorithms available in the software, was used to describe the Polycarbonate. The excellent fit of the calculated model to the measured spectra, achieved by the analysis software can be seen overlaid on the measured data. The model used was saved in the software database for use in analysis of subsequent measurements on this sample.




Figure 1. T and R for sample 1-Polycarbonate. Orientation 0°, s-polarisation.




Figure 2. n & k for sample 1-Polycarbonate. Orientation 0°, s-polarisation.




Figure 3. T and R for sample 1-Polycarbonate. Orientation 0°, p-polarisation.

Without rotating the sample, the same measurements were made but this time the incident light was selected at p-polarisation. The resulting spectrum (Fig. 3) and complex refractive indices (Fig. 4) are presented. A slight difference in refractive index between the p- and s-polarisation can be seen but it is not significant.




Figure 4. n&k for sample 1-Polycarbonate. Orientation 0°, p-polarisation.

The sample was then rotated by 90° and the measurements repeated. T and R for p-polarised incidence is shown in Figure 5. below. Analysis using the stored database model provided the refractive indices shown in Figure 6. Some difference is seen between this and the 0° sample orientation but again it is not significant.




Figure 5. T and R for sample 1-Polycarbonate. Orientation 90°, p-polarisation.




Figure 6. n&k for sample 1-Polycarbonate. Orientation 90°, p-polarisation.

The real changes are seen for the s-polarisation measurements at 90°, see Figure 7. Analysis of the spectrum, Figure 8, revealed a much higher refractive index than for all the other measurements. The degree of bi-refringence was calculated and is shown in Figure 9.




Figure 7. T and R for sample 1-Polycarbonate. Orientation 90°, s-polarisation.




Figure 8. n&k for sample 1-Polycarbonate. Orientation 90°, s-polarisation.




Figure 9. Difference between p- and s-polarisation measurements for the two measured orientations of the samples.

It is clear then that the nkd provides an ideal solution for measuring transparent and bi-refringent polymer films. In addition to providing precise transmittance and reflectance data the nkd can also be used to determine the refractive index, absorption and amount of bi-refringence in the material, quickly and accurately.